

STAR

Pediatric aerosol therapy system

medel[®]

medel[®] Baby



For babies



Easy-to-use

Medel **STAR** is a piston type compressor nebulizer system. The silent operation and the convenience of use make it an ideal device for home use.

Medel Star comes with a very attractive design, to ensure a pleasant therapy time for parents and children. It is equipped with a complete range of accessories: nebulizer, mouthpiece, adult mask, pediatric mask, nose-piece and air tube.

A funny colouring book with coloured pencils are included!

Class IIa Medical Device.



COLOURING BOOK WITH PENCILS

FAST

2ml in 5 minutes
only!

52
dBA

Silent: suitable for babies



Medeljet Basic



TECHNICAL SPECIFICATIONS

Flow	5,3 l/min
Aerosol output	0,171 ml
Aerosol output rate	0,082 ml/min
Nebulization rate	0,4 ml/min
Particle size (MMAD)	3,16 μ m
Fine particle fraction (FPF)	71,20%
Maximum compressor pressure	2,2 bar
Noise level	52 dBA
Minimum filling volume	2 ml
Maximum filling volume	8 ml
Operating mode	20 min. on / 40 min. off
Weight	1,5 Kg
Size	130 x 130 x 117mm



Accessories bag included!

95141

