

Fingertip Pulse Oximeter

Model: SB100




- %SpO2 & Pulse Rate - Big LED screen display for easy reading
- Shielded Design blocks Ambient Light
- Biocompatibility & Anti-Allergic Design
- One-touch ON/OFF key
- Pulse Strength Indicator
- Cord-attached - Economical & portable health management tool



%SpO2 & Pulse Rate



Shielded Design blocks Ambient Light



Biocompatibility & Anti-Allergic Design



Pulse Strength Indication



String-Attached



Large Font Size

Specification

Measurement Range | SpO2: 35%~99%; PR: 30~250 bpm

Precision | SpO2: 70%~99%; ±2%; 35%~69% (unspecified); Pulse rate: 30~250bpm; ±3

SpO2 Accessory | Adult Probe Included, Pediatric / Neonatal Probe Optional

Power Supply | 2 x "AAA" Alkaline

Pulse Oximeter

- Ideal for oxygen therapy assessment
- Ideal for sleep apnea screening
- Ideal for health management
- Ideal for spot-check & continuous monitoring

Model	Qty per carton	Carton volume
SB100	40 pcs	0.028 cbm/ctn

- Storage temperature: -25°C - 70°C (-13°F - 158°F)
- Relative humidity: 15% - 90% (no condensing)
- Atmospheric Pressure: 700 to 1060 hPa